



Mindful

Mental Health Support Chatbot

User Guide

How to use the application, end to end

Version 1.0 - June 2026

1. About Mindful

Mindful is a cross-platform mobile app that offers first-line emotional wellness support through a friendly chatbot. You can talk to it by typing or by speaking, and it can reply in text and read its replies aloud in a natural voice. Mindful also provides self-help resources and gentle coping tips.

Important: Mindful provides supportive, first-line help only. It is not a therapist and does not give medical diagnoses. If you are in crisis, please use the helpline information the app shows and reach out to a real person.

Key features

- Empathetic AI chat that understands how you feel.
- Voice input: speak your message and have it transcribed.
- Voice output: the assistant reads replies aloud (male or female voice).
- Emotion detection shown on each reply (e.g. stressed, sad, positive).
- Automatic crisis support with helpline information.
- Browse self-help resources by topic.
- Your conversations are stored only on your device (privacy-first).
- Admin panel to manage users and review chatbot logs.

2. Getting Started

What you need

- A smartphone (Android or iPhone).
- An internet connection for the chat and voice features.

Opening the app

Install Mindful on your device and tap its icon to open it. When you open Mindful for the first time, you will see a welcome screen with the logo, followed by the Login screen.

3. Creating an Account

1. On the Login screen, tap 'Register' at the bottom.
2. Enter your Full Name, Email, a Password (at least 6 characters), and confirm the password.
3. Tap the 'Register' button.
4. You are signed in automatically and taken to the Home screen.

Your account information is kept private and secure on your device.

4. Logging In

1. Enter the Email and Password you registered with.
2. Tap 'Login'.
3. If you forget your password, tap 'Forgot Password?'.

Mindful remembers you, so the next time you open the app it takes you straight in without logging in again - until you log out.

5. Home / Dashboard

The Home tab greets you by name and shows a 'Tip of the day'. Below it are quick-access cards:

- Talk to Mindful - opens the chat.
- Browse Resources - opens the self-help resources.
- Your Profile - opens your account and settings.

At the bottom of the screen is the navigation bar with four tabs: Home, Chat, Resources, and Profile.

6. Talking to Mindful (Chat)

Typing a message

1. Open the Chat tab.
2. Type how you feel in the message box at the bottom.
3. Tap the send button. Mindful replies with a supportive message.

Speaking instead of typing (voice input)

1. Tap the microphone icon. Recording starts and a moving waveform appears.
2. Speak your message.
3. Tap the check (tick) button when you are done.
4. Your speech is transcribed into the message box - review or edit it, then tap send. (Tap the X to cancel a recording.)

Hearing replies (voice output)

- Tap the speaker icon on any of Mindful's messages to hear it read aloud.
- While it speaks, an animated indicator shows on that message. Tap it again to stop.
- If you sent your message by voice, the reply is read aloud automatically.

Choosing the voice

- Tap the voice icon at the top of the Chat screen.
- Choose 'Female voice' or 'Male voice'. Your choice is saved.

Starting a new conversation

Tap the '+' button at the top right of the Chat screen to begin a fresh conversation. Your previous conversation stays saved, and the app reopens your most recent chat the next time you return.

7. Safety and Crisis Support

Mindful watches for messages that suggest you may be in serious distress (for example, thoughts of self-harm). If it detects this, a red support banner appears at the top of the chat with helpline contacts.

This check runs on your device and works even if the internet is down. Please reach out to the listed helplines or a trusted person if you see this banner.

8. Browse Resources

1. Open the Resources tab.
2. Use the topic chips at the top (All, Anxiety, Stress, Depression, Sleep, Mindfulness) to filter the list.
3. Tap any resource to read a short summary.

9. Profile and Logout

- Open the Profile tab to see your name, email, and account type.
- Tap 'Log out' to sign out. You will return to the Login screen.

10. Admin Panel (Administrators)

Administrators have an extra 'Admin Panel' option on the Profile screen. An administrator account is provided so an administrator can sign in straight away.

Administrator login

Email: admin@mindful.app

Password: Admin@1234

Using the Admin Panel

1. Log in with the administrator account above.
2. Go to the Profile tab and tap 'Admin Panel'.
3. Users tab: see everyone who has registered; remove a user with the delete icon.
4. Chatbot Logs tab: review recent chatbot messages, including the detected emotion and any crisis flags.

11. Privacy

Mindful is designed with your privacy in mind. Your conversations and account are kept only on your device and are not uploaded to any online account. Your messages are used only to generate helpful replies.

12. Troubleshooting / FAQ

The assistant does not reply

Check your internet connection and try again.

No sound when I tap the speaker

Make sure your phone is not on silent and the volume is up. Voice output needs an internet connection.

The microphone does not record

Allow microphone access when your phone asks for it, then try recording again.

I forgot my password

Use the 'Forgot Password?' option on the login screen, or contact your administrator for assistance.